




# IMS 1300-HR<sup>3</sup> KLEORA

## The Specialized Ion Microprobe for World-leading Research in Geochronology



- Best SIMS analytical performance
- Highest precision isotopic measurements
- Spot analysis, imaging, depth profiling
- Ease-of use & high throughput


In order to meet a growing demand from geochronologists, CAMECA introduces KLEORA, a new SIMS instrument fully optimized for advanced mineral dating.

Based on the new generation ultra-high sensitivity large geometry IMS 1300-HR<sup>3</sup> ion microprobe, KLEORA provides benchmark sensitivity for in-situ U-Th-Pb isotopic analyses in a high throughput, easy-to-use platform.

Both the IMS 1300-HR<sup>3</sup> and KLEORA inherit the instrumental excellence of CAMECA's ultra high sensitivity ion microprobes already adopted by top-ranked geoscience labs.

#### Key features inherited from the IMS 1280-HR and previous models:

- Large geometry design to provide high transmission at high mass resolution, required for analyzing low Pb concentration while eliminating interfering molecular ions.
- Oxygen flooding technique for improved Pb sensitivity (up to a factor of 7) and highly reproducible analytical conditions
- Superior ion imaging capabilities for trace element mapping in zircon grains (inhomogeneity, zonation)
- Enhanced magnet control system for high reproducibility at high mass resolution
- Remote operation, full automation, powerful geochronology dedicated software




#### Ion imaging for problematic zircon dating

Top: Pb and Ti scanning ion images.  
Bottom:  $^{207}\text{Pb}/^{206}\text{Pb}$  ratio image revealing a surprising  $\mu\text{m}$ -scale patchy Pb ratio distribution (Pb/Pb ages for ellipse areas in Ma).  
Data from M.A. Kusiak et al., *Geology* (2013)

## KLEORA benefits from the breakthrough instrumental advances of the IMS 1300-HR<sup>3</sup>

- **High brightness RF-plasma oxygen ion source** with greatly enhanced beam density and current stability, dramatically improving spatial resolution, data reproducibility and throughput **new source**
- **Automated sample loading system with motorized sample height (Z) adjustment** significantly increasing analysis precision, ease-of-use and productivity **auto storage** **auto-z**
- **UV-light microscope** for enhanced optical image resolution, together with dedicated software for easy sample navigation (developed by University of Wisconsin, USA) **uv-light**

With KLEORA, CAMECA offers a high throughput, easy-to-use ion microprobe, uniquely combining High Reproducibility with High spatial Resolution and High mass Resolution, and perfectly tailored to the needs of geochronologists!



#### U-Pb dating in 91500 zircon using RF-plasma source.

Excellent age precision, better than 0.3% ( $2\sigma$ ,  $n = 52$  spots) at high spatial resolution ( $\text{O}_2^-$  projected primary beam,  $I_p = 20\text{nA}$ , beam size  $< 10\ \mu\text{m}$ ).

High throughput (5 analysis per hour) with no compromise on data reproducibility.

IMS 1300-HR<sup>3</sup>: High Reproducibility • High spatial Resolution • High mass Resolution


## KLEORA delivers:

- In situ SIMS analysis with excellent sensitivity at high spatial resolution
- Superior analytical performance: isotopic ratio precision down to permil level
- Versatility: spot analysis, ion imaging, depth profiling...
- Ease-of-use & high throughput

## KLEORA is the ultimate solution for U-Pb dating due to:


- **Benchmark precision:** better than 0.3% ( $2\sigma$ ) for age precision in 91500 standard zircon.
  - Excellent data statistics (reduced uncertainty) due to high density primary beam **new source**
  - Benchmark long-term repeatability due to highly stable primary current **new source**
  - Highly reproducible conditions due to oxygen flooding and automated sample height adjustment **auto-z**
- **High productivity:** typically 5 analyses/hour for U-Pb spot analyses.
  - Short analysis time while keeping high precision and lateral resolution **new source**
  - Reduced maintenance (ion source and primary beam apertures with longer lifetime) **new source**
  - Multiple sample mounts analyzed in automated, unattended mode **auto storage**
- **Small spot size with high beam density** for excellent lateral resolution: tens of micron beam size with  $> 100 \text{ mA/cm}^2$  density, down to  $0.5\mu\text{m}$  with  $> 10 \text{ mA/cm}^2$  **new source**
- **Outstanding ion imaging capabilities:** sub-micron lateral resolution for microscope mode, and better than  $0.5\mu\text{m}$  for scanning ion imaging **new source**
- **Optimized depth profiling** with high sputtering rate, excellent sensitivity and dynamic range
- **Ease-of-use:**
  - computer-controlled sample exchange and sample height adjustment **auto storage** **auto-z**
  - easy analysis area selection **uv-light**
  - remote control, multi-user operation
- **Dedicated software:**
  - Geochronology data reduction module for U-Pb age determination
  - Winimage and Wincurve for powerful imaging and depth profiling data processing
- **Fast installation:** pre-tuned at our factory, the system is typically commissioned in less than six weeks
- **Low running cost** due to limited consumables and reduced maintenance of the new RF-plasma source
- **Standards:** includes geochronology standard samples.

High spatial resolution, in situ geochronology (U-Pb, U-Th), either grain by grain, or by subgrain domains, revealing complexities that cannot be resolved by bulk analysis methods.



Magmatic dating and isotope study of individual crystals from the Yellowstone supervolcano:  
CL image of zircon dated for U-Th ages (blue) and U-Pb ages (yellow).  
Data from K. E. Watts et al., Contrib Mineral Petrol (2012)

Unique adjustable, small spot size capabilities using Gaussian primary ion beam illumination instead of conventional projected mode.



U-Pb dating at a scale  $< 5\mu\text{m}$  on zircon standards demonstrating close to 1% age precision and accuracy. Gaussian primary illumination mode,  $O_2$  spot size  $< 5\mu\text{m}$ .

Data from Y. Liu et al., JAAS (2011)

### Standard configuration

	IMS 1300-HR <sup>3</sup>	KLEORA
RF-plasma O source <b>New</b>	•	•
Cs microbeam source	•	
Auto-Z motion <b>New</b>	•	•
Auto storage chamber <b>New</b>	•	•
High resolution UV- light microscope <b>New</b>	•	•
Oxygen flooding	•	•
Normal-incidence Electron Gun	•	
Multicollector system	•	
NMR control	•	
Geochronology software	•	•

**KLEORA can be upgraded to a full IMS 1300-HR<sup>3</sup> model at any time.**




Based on the new IMS 1300-HR<sup>3</sup> ion microprobe, **KLEORA** delivers superior SIMS analytical performance in a high throughput, easy-to-use platform, covering the entire range of geochronological applications.

## Powerful U-Pb dating data processing software

KLEORA includes a specific geochronology data reduction software running on MS Excel™ platform, developed at NORDSIM by Dr. Martin Whitehouse.

It offers outstanding features for accurate and highly productive data processing:

- Direct import of SIMS files in Excel
- Fast and easy processing of large volumes of data
- Multiple types of calibration including Pb/U vs. UO<sub>2</sub>/U using power law parameters
- Common lead correction methods
- Age table display
- Export for Concordia diagram plotting.



From precise age determination of zircon to dating of other U-rich minerals such as apatite, rutile or baddeleyite, **KLEORA** provides new types of innovative dating methods. Refer to our Geochronology application note for more information.

**CAMECA** is the world premier provider of microanalytical instrumentation. We deliver cutting-edge science and metrology solutions, and offer our customers unparalleled support and maintenance service through the comprehensive AMECARE program.

### Corporate Headquarters:

29 Quai des Grésillons • 92622 Gennevilliers Cedex • FRANCE • Phone: +33 1 43 34 62 00 • [cameca.info@ametek.com](mailto:cameca.info@ametek.com)

### CAMECA Subsidiaries or International AMETEK Offices:

AMETEK do Brasil CAMECA Business Unit • [cameca@ametek.com.br](mailto:cameca@ametek.com.br)  
 AMETEK China CAMECA Business Unit • [cameca-china.sales@ametek.com.cn](mailto:cameca-china.sales@ametek.com.cn)  
 CAMECA Germany • [sales.cameca@ametek.de](mailto:sales.cameca@ametek.de)  
 AMETEK India CAMECA Business Unit • [cameca-india@ametek.com](mailto:cameca-india@ametek.com)  
 AMETEK Japan CAMECA Business Unit • [cameca-japan.info@ametek.co.jp](mailto:cameca-japan.info@ametek.co.jp)  
 AMETEK Korea CAMECA Business Unit • [cameca-korea.info@ametek.co.kr](mailto:cameca-korea.info@ametek.co.kr)  
 CAMECA Taiwan Corp. Ltd. • [sales.cameca@ametek.com.tw](mailto:sales.cameca@ametek.com.tw)  
 CAMECA Instruments, Inc. USA • [cameca.us-sales@ametek.com](mailto:cameca.us-sales@ametek.com)

